



PELI
BIO THERMAL™

EMT



EMT, FEATURING GOLDEN HOUR™ TECHNOLOGY

Based on our award winning Golden Hour™ technology, the EMT was developed to be used in emergency situations by first responders, emergency medical services and paramedics. EMT can be used to safely transport blood to an accident site in a critical amount of time.

The very same patented and award-winning technology provides superior thermal protection for high-value, temperature-sensitive payloads in our renowned Military products; Golden Hour™ Mobile and Golden Hour™ Medic containers.

WHY CHOOSE EMT SYSTEMS?

- Ease of use - modular design for quick assembly and efficient inventory storage
- Low profile and life-saving
- Nimble and tactical for delivery of whole blood to the point of injury within a crucial window of time
- Quality, durable construction
- Designed and tested to last in the harshest conditions

EMT



DURATION

48 - 72 hours



TEMPERATURE

SERIES 4

2°C to 8°C



PAYLOAD SIZE

2L




EMT PRODUCT SPECIFICATION

MOBILE	Payload Size	Litres	Outer Material	Exterior Dimensions WxDxH		Interior Dimensions WxDxH		Tare Weight	
				mm	inches	mm	inches	kg	lbs
				Series 4 EMT	2L	Standard Bag	235 x 210 x 203	9.25 x 8.25 x 8	157 x 127 x 108

*Fully loaded container with blood product payload

ADDITIONAL DATA



FOR MORE INFORMATION ABOUT PELI BIOTHERMAL PRODUCTS [CLICK HERE](#)

FOR MORE INFORMATION ABOUT PELI BIOTHERMAL PRODUCTS, PLEASE CONTACT:

EMEA OFFICE - INFOEMEA@PELIBIOTHERMAL.COM T: +44 (0)1525 243 770

AMERICAS OFFICE - INFOAMERICAS@PELIBIOTHERMAL.COM T: +1(763) 412-4800 TOLL FREE US: T: +(877) 537-9800

ASIA OFFICE - INFOASIA@PELIBIOTHERMAL.COM T: +65 6681 0095

PELIBIOTHERMAL.COM

ISO17025, ISO9001-2008 Accredited

Golden Hour™, TIC™, Crêdo™, Crêdo Cube™, Crêdo ProEnvision™, Crêdo Xtreme™, Crêdo Duracube™, Crêdo ProMed™, Sherpa Systems™, CoolGuard™ and CoolPall™ are registered trademarks of Pelican Products Inc., its subsidiaries and / or affiliates



©2022 Peli BioThermal Ltd all rights reserved